

NEVADA DEPARTMENT OF PUBLIC SAFETY

BREATH TEST OPERATOR INTOXILYZER 8000 CHECKLIST

INSTRUMENT SERIAL # _____

AGENCY: _____

CASE #: _____

SUBJECT _____

DATE: _____

OPERATOR: _____

CERTIFICATION #: _____

If the instrument displays – “STANDBY MODE” – press the green START TEST button.
(Note: Instrument displays a Ready Mode count down and goes through diagnostics before it is ready to begin.)

1. Check the subject's mouth for foreign objects, including dental work and studs or other jewelry, and have the subject remove any readily removable objects. Make a note of any permanent objects or foreign objects not removed.
2. **NOTE THE TIME:** _____, Begin a deprivation period of at least 15 minutes. Once started, observe the subject and do not allow him/her to eat, drink, smoke or put any foreign object in his/her mouth. If the subject vomits or appears to regurgitate, and you wish to continue with breath testing, you must begin a new 15 minute deprivation period.

Note any observations regarding objects in the mouth or the deprivation period (or any other relevant comments) here:

It is recommended that at this point you ensure adequate paper is in internal printer feeder or external printer paper tray.

3. Verify the instrument displays the “READY MODE” message and push the green START TEST button to start.
4. When the display requests “SWIPE/SCAN OPERATOR ID OR PRESS ENTER”, place the barcode of the Operator Card under the red LED scanner. If accepted, press ENTER to review scanned information. If the barcode is not accepted or not available, press ENTER and answer the subsequent questions, followed by ENTER.
5. When the display requests “PLEASE SWIPE/SCAN DL OR PRESS ENTER”, place the barcode of the subject's driver's license under red LED scanner or swipe the magnetic stripe using the reader on the top of the instrument. If accepted, press ENTER to review the uploaded subject information. Then answer the additional test data entry questions, followed by ENTER. If the barcode is not accepted, press ENTER and answer the subsequent test data entry questions, followed by ENTER.
6. The display will request, “OBSERVATION TIME?” Enter the time that the deprivation period began (time listed on Step #2) followed by ENTER. The instrument will automatically run an air blank, diagnostic test, air blank and a standard test. If the standard test is out of range, the instrument will not allow the completion of the subject's breath test.
7. When the display requests “PLEASE BLOW UNTIL TONE STOPS / R” attach a clean mouthpiece and request the subject blow with a long, continuous breath into the breath tube until the tone stops. If the subject is not willing to provide a sample, you may press the “R” key. The instrument will only accept this command when “PLEASE BLOW UNTIL TONE STOPS / R” is on the display. After accepting the sample the display will then show, “PLEASE WAIT” as it counts down the two minute wait period before the second subject test sample.
8. When the display requests “PLEASE BLOW UNTIL TONE STOPS / R” attach a clean mouthpiece and request the subject blow with a long, continuous breath into the breath tube until the tone stops. If the subject is not willing to provide a sample, you may press the “R” key. The instrument will only accept this command when “PLEASE BLOW UNTIL TONE STOPS / R” is on the display.
9. (Conditional Step) If the measured results of the two samples taken above do not agree within 0.020, the instrument will automatically request a third sample be given after a two minute wait period. When requested, have the subject deliver a third sample as described in steps 8 and 9 above.
10. Upon completion of the testing sequence, the instrument will automatically print out the test results. **REMOVE THE TEST PRINTOUT and SIGN AND DATE IT. CORRECT ANY MISINFORMATION ON THE EVIDENCE TEST PRINTOUT IF NECESSARY and INITIAL THE CHANGES.** RECORD required information in the D.U.I. LOGBOOK. For multiple copies, press F2.

I HAVE FOLLOWED THE PROCEDURES OUTLINED ABOVE.

OPERATOR'S SIGNATURE